



Anti-CD79A polyclonal antibody [Texas Red®] (CABT-BL1006)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Human IgA alpha heavy chain
Source/Host	Goat
Species Reactivity	Human
Purification	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Human antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities, pepsin digestion and chromatographic separation
Conjugate	Texas Red
Applications	IF
Format	Lyophilized
Buffer	0.01 M Sodium Phosphate, 0.25 M Sodium Chloride, pH 7.2
Preservative	None
Storage	Store vial at 4°C prior to restoration. For extended storage aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature.

BACKGROUND

Introduction	ELF3 (E74-like factor 3 (ets domain transcription factor, epithelial-specific)) is a protein-coding gene. Diseases associated with ELF3 include retroperitoneal sarcoma, and pandas. GO annotations related to this gene include transcription coactivator activity and sequence-specific
---------------------	--

DNA binding transcription factor activity. An important paralog of this gene is ELK4.

GENE INFORMATION

Entrez Gene ID	973
----------------	---------------------

Protein Refseq	NP_001774
----------------	---------------------------

UniProt ID	P11912
------------	------------------------
