



# Anti-CDC11 (C-terminal) polyclonal antibody (CABT-BL1010)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide derived from the C terminal domain of CDC11
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	<i>S. cerevisiae</i>
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, WB
<b>Cellular Localization</b>	Membrane. Bud neck.Note: Present at the bud neck during cell division.
<b>Format</b>	Liquid
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">853539</a>
<b>Protein Refseq</b>	<a href="#">NP_012610</a>

UniProt ID

[P32458](#)

---