



Anti-CDC3 (N-terminal) polyclonal antibody (CABT-BL1013)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide derived from N terminal domain of CDC3 (<i>Saccharomyces cerevisiae</i>).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Cell Membrane
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	851024
Protein Refseq	NP_013418

UniProt ID

[P32457](#)
