



Anti-DUXBL (aa 241-290) polyclonal antibody (CABT-BL1326)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: PTTTPGEQGF HSQPPLQLLT QNRGHNPRESGGLAVPRLED CTQVPAVNQH, corresponding to amino acids 241-290 of Mouse Duxbl
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	278672
Protein Refseq	NP_899245

