



Anti-ELA1 (aa 115-214) polyclonal antibody (CABT-BL1362)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Fusion protein: VKDIRNNKYRIPYRMLYSKYQQEVEKKQEESAERLRLEMQ KLQQEREKKQTIVVDHTVYFKKRNTKKTTTRLHNEPHSQLY MKSLKDHESRLKHFKDGGFN, corresponding to amino acids 115/214 of <i>S. cerevisiae</i> Elongin A
Isotype	IgG
Source/Host	Mouse
Species Reactivity	<i>S. cerevisiae</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	50% Glycerol
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Acts upon elastin.
---------------------	--------------------

GENE INFORMATION

Entrez Gene ID [855491](#)

Protein Refseq [NP_014169](#)

UniProt ID [P53861](#)
