



Anti-ELFN2 (aa 254-303) polyclonal antibody (CABT-BL1364)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 254-303 (ARPVSHPTPY STDAQREPDE NSGFNPDEIL SVEPPASSTT DASAGPAIKL) of human LRRC62.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Membrane; Single pass membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [114794](#)

Protein Refseq [NP_443138](#)

UniProt ID [Q5R3F8](#)
