



# Anti-EML3 (aa 180-229) polyclonal antibody (CABT-BL1369)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                                                                 |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | A synthetic peptide corresponding to a region within the N terminal amino acids 180-229 (LLVRSGSTES RGGKDPLSSP GGPGSRRSNY NLEGISVKMF LRGRPITMYI) of human EML3. |
| <b>Isotype</b>               | IgG                                                                                                                                                             |
| <b>Source/Host</b>           | Rabbit                                                                                                                                                          |
| <b>Species Reactivity</b>    | Human                                                                                                                                                           |
| <b>Purification</b>          | Immunogen affinity purified                                                                                                                                     |
| <b>Conjugate</b>             | Unconjugated                                                                                                                                                    |
| <b>Applications</b>          | WB                                                                                                                                                              |
| <b>Cellular Localization</b> | Cytoplasmic                                                                                                                                                     |
| <b>Format</b>                | Liquid                                                                                                                                                          |
| <b>Buffer</b>                | 2% Sucrose, PBS                                                                                                                                                 |
| <b>Preservative</b>          | None                                                                                                                                                            |
| <b>Storage</b>               | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.                                                                  |

## GENE INFORMATION

Entrez Gene ID [256364](#)

Protein Refseq [NP\\_694997](#)

---

UniProt ID [Q32P44](#)

---