



Anti-EPHX3 polyclonal antibody (CABT-BL1381)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	The antiserum was produced against synthesized peptide derived from internal of human ABHD9.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C/1 year

BACKGROUND

Introduction	EPHX3 (Epoxide Hydrolase 3) is a Protein Coding gene. GO annotations related to this gene include hydrolase activity. An important paralog of this gene is EPHX4.
---------------------	---

GENE INFORMATION

Entrez Gene ID [79852](#)

Protein Refseq [NP_001136358](#)

UniProt ID [A0A024R7F3](#)
