



# Anti-ERCC4 (aa 855-905) polyclonal antibody (CABT-BL1387)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region between residues 855 and 905 of human XPF (AAB07689.1)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Cellular Localization	Nucleus.
Format	Liquid
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

## BACKGROUND

Introduction	The protein encoded by this gene forms a complex with ERCC1 and is involved in the 5' incision made during nucleotide excision repair. This complex is a structure specific DNA repair endonuclease that interacts with EME1. Defects in this gene are a cause of xeroderma
--------------	---

## GENE INFORMATION

Entrez Gene ID	<a href="#">2072</a>
----------------	----------------------

---

Protein Refseq	<a href="#">NP_005227</a>
----------------	---------------------------

---

UniProt ID	<a href="#">Q92889</a>
------------	------------------------

---