



Anti-ergo-1 (aa 1103-1121) polyclonal antibody (CABT-BL1393)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: CEVNKDMNVNEKLEGMTFV, corresponding to C terminal amino acids 1103-1121 of <i>Caenorhabditis elegans</i> ergo1. The peptide was amidated and conjugated to a carrier protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>Caenorhabditis elegans</i>
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	178602
Protein Refseq	NP_503362
UniProt ID	O61931