



Anti-ESYT3 (aa 72-121) polyclonal antibody (CABT-BL1398)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: WRRNRRGKLG RLAAAFEFLD NEREFISREL RGQHLPAAWIH FPDVERVEWA, corresponding to N terminal amino acids 72-121 of Human FAM62C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell membrane; Single-pass membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	83850
Protein Refseq	NP_114119

UniProt ID

[A0FGR9](#)
