



Anti-EVE (aa 1-50) polyclonal antibody (CABT-BL1402)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MHGYRTYNME SHHAHHDASP VDQKPLVVDL LATQYGKPQT PPPSPNECLS) of Human EVE.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	The eve gene contains a homeo box and hence is related structurally to the pair-rule gene fushi tarazu and to homeotic selector genes.
---------------------	--

GENE INFORMATION

Entrez Gene ID [36039](#)

Protein Refseq [NP_523670](#)

UniProt ID [F9W320](#)
