



Anti-EVI2A (full length) polyclonal antibody (CABT-BL1403)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein, corresponding to amino acids 1-237 of Human EVI2A
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Cell Membrane
Format	Liquid
Concentration	Lot specific
Buffer	1X PBS, pH 7.2
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [2123](#)

Protein Refseq [NP_001003927](#)

UniProt ID [P22794](#)
