



Anti-EXD (aa 288-327) polyclonal antibody (CABT-BL1404)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within a synthetic peptide: WFGNKRIRYK KNIGKAQEEA NLYAAKKAAG ASPYSMAGPP SGTTPMMSP, corresponding to C terminal amino acids 288-327 of Fruit fly (Drosophila melanogaster) Homeobox protein extradenticle
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Concentration	Lot specific
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [32567](#)

Protein Refseq [NP_523360](#)

UniProt ID [M9PHV1](#)
