



Anti-EXPH5 (internal region) polyclonal antibody (CABT-BL1406)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 1907-1956 (QGRLWKPSFLKNPGFLKDDLNRNPPNPSESLSSNSPSSQV PEDGLSPSEPL) of Human EXPH5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Concentration	Lot specific
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Exophilin 5, also known as EXPH5, is a human gene.

GENE INFORMATION

Entrez Gene ID [23086](#)

Protein Refseq [NP_055880](#)

UniProt ID [Q149M6](#)
