



Anti-FAM101A (aa 73-122) polyclonal antibody (CABT-BL1412)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within the internal region amino acids 73-122 (LEPRPRALRF RSTTIIFPKH ARSTFRTLH CSLGRPSRWF TASVQLQLCQ) of Human FAM101A (Accession number: EAW98447)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Concentration	Lot specific
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [144347](#)

Protein Refseq [NP_859060](#)

UniProt ID [Q6ZTI6](#)
