



Anti-FAM103A1 polyclonal antibody (CABT-BL1413)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human FAM103A1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Concentration	Lot specific
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	FAM103A1 (family with sequence similarity 103, member A1) is a protein-coding gene. GO annotations related to this gene include protein binding and RNA binding.
--------------	--

GENE INFORMATION

Entrez Gene ID [83640](#)

Protein Refseq [NP_113640.1](#)

UniProt ID [Q9BTL3](#)
