



Anti-FAM105A (aa 181-230) polyclonal antibody (CABT-BL1416)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 181-230 (IQQYSFGPEK YTGSNVFGKL RKYVELLKTQ WTEFNGIRDY HKRGSMCNTL) of Human FAM105A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	54491
Protein Refseq	NP_061891

UniProt ID

[Q9NUU6](#)
