



# Anti-FAM109A (aa 35-84) polyclonal antibody (CABT-BL1417)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within N terminal amino acids 35-84 (HRRWFVLRGN MLFYFEDAAS REPVGVILLE GCTVELVEAA EEFAFAVRFA) of Human FAM109A.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	This gene encodes a protein that localizes to the endosome and interacts with the enzyme, inositol polyphosphate 5-phosphatase OCRL-1. Alternate splicing results in multiple transcript variants.
---------------------	--

# GENE INFORMATION

Entrez Gene ID [144717](#)

---

Protein Refseq [NP\\_001171467](#)

---

UniProt ID [Q8N4B1](#)

---