



Rabbit Anti-Human FAM113B Polyclonal Antibody (CABT-BL1419)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human FAM113B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	PCED1B;PC-esterase domain containing 1B;FAM113B;PC-esterase domain-containing protein 1B;family with sequence similarity 113, member B
----------	--

GENE INFORMATION

Gene Name	PCED1B
Entrez Gene ID	91523
Protein Refseq	NP_001071178.2
UniProt ID	A0A024R115