



Anti-FAM116B polyclonal antibody (CABT-BL1420)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human FAM123A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	DENND6B (DENN/MADD domain containing 6B) is a protein-coding gene. GO annotations related to this gene include Rab guanyl-nucleotide exchange factor activity. An important paralog of this gene is DENND6A.
--------------	--

GENE INFORMATION

Entrez Gene ID [414918](#)

Protein Refseq [NP_001001794.3](#)

UniProt ID [Q8NEG7](#)
