



Anti-FAM118A (aa 215-264) polyclonal antibody (CABT-BL1421)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 215-264 (EVMEVLQONLY RTKSFLFVGC GETLRDQIFQ ALFLYSVPNK VDLEHYMLVL) of Human FAM118A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [55007](#)

Protein Refseq [NP_060381](#)

UniProt ID [A0A024R4V3](#)
