



Anti-FAM129B polyclonal antibody (CABT-BL1429)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within amino acids 683-733 of Human FAM129B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IP
Format	Liquid
Size	100 µl
Buffer	0.1% BSA, Tris buffered saline
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	64855
Protein Refseq	NP_073744

UniProt ID

[A0A024R872](#)
