



Anti-FAM129B polyclonal antibody (CABT-BL1430)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Peptide sequence around aa. 710~714 (Q-V-S-S-P) derived from Human Niban-like protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IF, WB
Buffer	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C/1 year

GENE INFORMATION

Entrez Gene ID	64855
Protein Refseq	NP_001030611
UniProt ID	A0A024R872