



Anti-FAM131C (aa 107-156) polyclonal antibody (CABT-BL1431)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 107-156 (RGRVAHLIEW KGWSAQPAGW ELSPAEDEHY CCLPDELREA RFAAGVAEQF) of Human FAM131C
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	348487
Protein Refseq	NP_872429

UniProt ID

[Q96AQ9](#)
