



Anti-FAM135B (aa 1259-1308) polyclonal antibody (CABT-BL1432)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 1259-1308 (TLVSTGLWLM QKLKKSGLL QLTFRDNADL RKCFLYQLSQ KTGLQYFKNV) of Human FAM135B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	51059
Protein Refseq	NP_056996

UniProt ID

[Q49AJ0](#)
