



Anti-FAM154A (aa 36-85) polyclonal antibody (CABT-BL1436)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 36-85 (EYTENYPFYH SYLPRESFKP RREYQKGSIP MEGLTTSRRD FGPHKVAPVK) of Human FAM154A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	158297
Protein Refseq	NP_714918

UniProt ID

[F6S232](#)
