



Anti-FAM164C (aa 110-159) polyclonal antibody (CABT-BL1439)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from within a region of N terminal amino acids 110-159 (SSGPQSWYPK ANNQDFIPFT KKRVGVDRAF PLKPMVHRKS CSTGEAGTDG) of Human C14ORF140
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	79696
Protein Refseq	NP_078919
UniProt ID	A0A024R6E6