



# Anti-FAM167A polyclonal antibody (CABT-BL1441)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	Synthetic peptide selected from the central region of Human FAM167A, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID	<a href="#">83648</a>
Protein Refseq	<a href="#">NP_444509</a>
UniProt ID	<a href="#">Q96KS9</a>