



Anti-FAM171A1 (internal region) polyclonal antibody (CABT-BL1445)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 468-517 (ARKSMEREGYESSGNDYRGSYNTVLSQPLFEKQDREGP ASTGSKLTIQE) of Human FAM171A1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Membrane; Single-pass type I membrane protein (Potential)
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [221061](#)

Protein Refseq [NP_001010924](#)

UniProt ID [B3KMX9](#)
