



Anti-FAM174B (aa 109-158) polyclonal antibody (CABT-BL1446)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: LLRVFRSGKR LKKTRKYDII TTPAERVEMA PLNEEDDEDE DSTVFDIKYR, corresponding to C terminal amino acids 109-158 of Human LOC400451
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	400451
Protein Refseq	NP_997329
UniProt ID	Q3ZCQ3