



Anti-FAM176A polyclonal antibody (CABT-BL1447)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide from human transmembrane protein 166 (TM166, TMEM166) conjugated to an immunogenic carrier protein has been used as the immunogen. The antigen is homologous in rat and mouse.
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC, WB
Format	Lyophilised
Preservative	None
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. When reconstituting, glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.
Ship	This item will be shipped to you at ambient temperature in a lyophilised form.

GENE INFORMATION

Entrez Gene ID	84141
Protein Refseq	NP_115557.1

UniProt ID

[Q9H8M9](#)
