



Anti-FAM190A (aa 787-816) polyclonal antibody (CABT-BL1455)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, between amino acids 787-816 from the C-terminal region of Human FAM190A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	401145
Protein Refseq	NP_001138537
UniProt ID	Q9C0I3