



Anti-FAM190A (aa 72-121) polyclonal antibody (CABT-BL1456)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide designed within residues: HKKGSEPKQE PTNQNLISIN GAQPGHSNMQ KLSLEEHIKT RGRHSVGFSS, corresponding to amino acids 72 - 121 of Human FAM190A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	401145
Protein Refseq	NP_001138537

UniProt ID

[Q9C0I3](#)
