



Anti-FAM36A (aa 71-118) polyclonal antibody (CABT-BL1471)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 71-118 (LGCWFHCRYN YAKQRIQERI AREEIKKKIL YEGTHLDPER KHNGSSSN) of Human FAM36A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Cellular Localization	Membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [116228](#)

Protein Refseq [NP_932342](#)

UniProt ID [B3KM21](#)
