



Anti-FAM40B (aa 20-50) polyclonal antibody (CABT-BL1473)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 20-50 of Human FA40B (isoform 2), conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. Aliquot and store at -20°C long-term.

GENE INFORMATION

Entrez Gene ID	57464
Protein Refseq	NP_001127808
UniProt ID	Q9ULQ0