



Anti-FAM47A (aa 504-553) polyclonal antibody (CABT-BL1474)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 504-553 (LRSEPPKTRR TSSLRSEPPK TRRTSSLGPE PPKTRRVSSL RPELPKSRRV) of Human FAM47A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	158724
Protein Refseq	NP_981953

UniProt ID

[Q5JRC9](#)
