



Anti-FAM50B (aa 35-84) polyclonal antibody (CABT-BL1475)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to a region within N terminal amino acids 35-84 (IAEETILKSQ VDKRFSAHYD AVEAELKSST VGLVTLNDMK ARQEALVRER) of Human FAM50B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Reconstitution	Add 50 µl of distilled water.
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	26240
Protein Refseq	NP_036267

UniProt ID

[A0A024QZY3](#)
