



# Anti-FAM59B (C-terminal) polyclonal antibody (CABT-BL1477)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A 19 amino acid synthetic peptide from near the carboxy terminus of Human FAM59B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA
Format	Liquid
Buffer	PBS
Preservative	0.02% Sodium Azide
Storage	Store at +4°C.

## BACKGROUND

Introduction	GAREML (GRB2 associated, regulator of MAPK1-like) is a protein-coding gene. An important paralog of this gene is GAREM.
--------------	---

## GENE INFORMATION

Entrez Gene ID [150946](#)

---

Protein Refseq [NP\\_001161713](#)

---

UniProt ID [Q75VX8](#)

---