



Anti-FAM71D (internal region) polyclonal antibody (CABT-BL1480)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 288-337 (TVVFENNDLIRAKQEEKEKLKNILKPGCLQDTKSKSELK ESSKHVTISNI) of Human FAM71D
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	161142
Protein Refseq	NP_775797
UniProt ID	Q8N9W8