



# Anti-FAM76B (internal region) polyclonal antibody (CABT-BL1483)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within internal sequence amino acids 108-157 (QQCAFDRKEE GRRKVDGKLL CWLCTLSYKR VLQKTKEQRK SLGSSHSNSS) of Human FAM76B
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, WB
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction** Protein FAM76B is a protein that in humans is encoded by the FAM76B gene.

# GENE INFORMATION

Entrez Gene ID [143684](#)

---

Protein Refseq [NP\\_653265](#)

---

UniProt ID [Q5HYJ3](#)

---