



# Anti-FAM81A (aa 72-121) polyclonal antibody (CABT-BL1484)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                           |                                                                                                                                                               |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>          | Synthetic peptide corresponding to a region within the N terminal amino acids 72-121(DRLARLFLEE HIRNITAIVK QLNRDIEVLQ EQIRARDNIS YGTNSALKTL) of HUMAN FAM81A. |
| <b>Isotype</b>            | IgG                                                                                                                                                           |
| <b>Source/Host</b>        | Rabbit                                                                                                                                                        |
| <b>Species Reactivity</b> | Human                                                                                                                                                         |
| <b>Purification</b>       | Immunogen affinity purified                                                                                                                                   |
| <b>Conjugate</b>          | Unconjugated                                                                                                                                                  |
| <b>Applications</b>       | WB                                                                                                                                                            |
| <b>Reconstitution</b>     | Add 50 µl of distilled water.                                                                                                                                 |
| <b>Format</b>             | Lyophilised                                                                                                                                                   |
| <b>Buffer</b>             | 2% Sucrose, PBS                                                                                                                                               |
| <b>Preservative</b>       | None                                                                                                                                                          |
| <b>Storage</b>            | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.                                                       |

## GENE INFORMATION

Entrez Gene ID [145773](#)

Protein Refseq [NP\\_689663](#)

---

UniProt ID [Q8TBF8](#)

---