



Anti-FAM81B (aa 2-51) polyclonal antibody (CABT-BL1485)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the N terminal amino acids 2-51 (QLQLGLTLAS SEKRKKSQRL FFKNIKSTKN KAGKASIMSS DTNVNKSASP) of Human FAM81B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	153643
Protein Refseq	NP_689761

UniProt ID

[Q96LP2](#)
