



# Anti-FAM82A1 polyclonal antibody (CABT-BL1486)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	This antibody is purified through a protein A column, followed by peptide affinity purification.
Conjugate	Unconjugated
Applications	WB, ELISA
Buffer	PBS with 0.09% (W/V) sodium azide
Preservative	0.09% Sodium Azide
Storage	RMD2 Antibody (Center) can be refrigerated at 2-8°C for up to 6 months. For long term storage, place the at -2°C.

## GENE INFORMATION

Entrez Gene ID	<a href="#">151393</a>
Protein Refseq	<a href="#">NP_653314.2</a>
UniProt ID	<a href="#">Q96LZ7</a>