



Anti-FAM83E (aa 359-408) polyclonal antibody (CABT-BL1487)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to internal sequence amino acids 359-408 (RARTPSGPPA RPSRSMWDLS RLSQLSGSSD GDNELKKS WG SKDTPAKALM) of Human FAM83E
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	54854
Protein Refseq	NP_060178
UniProt ID	Q2M2I3