



Anti-FAM83F (aa 143-192) polyclonal antibody (CABT-BL1488)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 143-192 (KDEKAPHLKQ VVRQMIQQAQ KVIADVMDLF TDGDIFQDIV DAACKRRVPV) of Human FAM83F.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	113828
Protein Refseq	NP_612444

UniProt ID

[Q8NEG4](#)
