



Anti-FAM91A1 (aa 300-350) polyclonal antibody (CABT-BL1489)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region between residue 300 and 350 of human FAM91A1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Format	Liquid
Size	100 µl
Buffer	8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8
Preservative	0.09% Sodium Azide
Storage	Store at +4°C.

GENE INFORMATION

Entrez Gene ID	157769
Protein Refseq	NP_659400

UniProt ID

[Q658Y4](#)
