



# Anti-FAT1 polyclonal antibody (CABT-BL1498)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Recombinant protein corresponding to amino acids of human FAT1.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Antigen affinity purification
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC
<b>Format</b>	Liquid
<b>Buffer</b>	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	This gene is an ortholog of the Drosophila fat gene, which encodes a tumor suppressor essential for controlling cell proliferation during Drosophila development. The gene product is a member of the cadherin superfamily, a group of integral membrane proteins characterized by the presence of cadherin-type repeats. In addition to containing 34 tandem cadherin-type repeats, the gene product has five epidermal growth factor (EGF)-like repeats and one laminin A-G domain. This gene is expressed at high levels in a number of fetal epithelia. Its product
---------------------	---

probably functions as an adhesion molecule and/or signaling receptor, and is likely to be important in developmental processes and cell communication. Transcript variants derived from alternative splicing and/or alternative promoter usage exist, but they have not been fully described. [provided by RefSeq, Jul 2008]

---

## GENE INFORMATION

Entrez Gene ID	<a href="#">2195</a>
----------------	----------------------

Protein Refseq	<a href="#">NP_005236.2</a>
----------------	-----------------------------

UniProt ID	<a href="#">Q14517</a>
------------	------------------------

---