



# Anti-FER1L5 (N-terminal) polyclonal antibody (CABT-BL1509)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide derived from the N terminal region of Human FER1L5.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC
<b>Cellular Localization</b>	Cell Membrane; single pass membrane protein
<b>Format</b>	Liquid
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">90342</a>
<b>Protein Refseq</b>	<a href="#">NP_001106853</a>

UniProt ID

[A0A096LNV2](#)

---