



Anti-FEZF1 (C-terminal) polyclonal antibody (CABT-BL1511)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide of the C terminal part of Mouse Fezf1 (CGPLPSGLHQGHQ)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Purified with antigen peptide
Conjugate	Unconjugated
Applications	IHC, WB
Format	Lyophilized product from 1 % BSA in PBS containing 0.05% NaN ₃ .
Preservative	0.05% Sodium Azide
Storage	Lyophilized product, 5 years at 2 - 8°C; Solution, 2 years at -20°C.

BACKGROUND

Introduction	Fasciculation and elongation protein zeta-1 is a protein that in humans is encoded by the FEZF1 gene.
--------------	---

GENE INFORMATION

Entrez Gene ID	73191
----------------	-----------------------

Protein Refseq [NP_082738](#)

UniProt ID [Q0VDQ9](#)
